| | |
|--|--|
|  2 - 3 | <i>MANȘETĂ DE PRESIUNE CU MANOMETRU OB</i> <i>INSTRUCȚIUNI DE UTILIZARE</i> |
|  4 - 5 | <i>OB PRESSURE INFUSION CUFF OB</i> <i>USER MANUAL</i> |
|  6 - 7 | <i>OB DRUCKINFUSIONSMANSCHETTE MIT MANOMETER</i> <i>BEDIENUNGSANLEITUNG</i> |

CE

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OB PRESSURE INFUSION CUFF

1. DESCRIPTION OF DEVICE

The OB pressure infusion cuff with manometer is a class I medical device as per European Directive MDD93/42/EEC. This device is not a measurement instrument. Its manometer only has a mere indicative purpose and allows you to know the pressure value applied on the infusion cuff.

The OB pressure infusion cuff with manometer is designed and manufactured to ensure a practical, safe and correct use. Please read this leaflet carefully before use. In case of doubts do not hesitate to contact the Oscar Boscarol srl.

2. TECHNICAL FEATURES

| | |
|----------------------------|--|
| Name | OB PRESSURE INFUSION CUFF WITH MANOMETER |
| Risk class | class I medical device - Directive 93/42/CEE |
| Dimensions | 33x17 cm |
| Pressure range | 0 ÷ 300 mmHg (millimeter of mercury) |
| Resolution | 2 mmHg |
| Precision | ±3 mmHg (temperature: 18°C ÷ 35°C) ±6 mmHg (temperature: 36°C ÷ 50°C) |
| Air lack | <4 mmHg/min |
| Use temperature | 0 ÷ 50 °C r.h. 10 ÷ 85 % |
| Storage temperature | -20 ÷ 70 °C r.h. 10 ÷ 85 % |

The data indicated are not binding for the manufacturer and may change without prior notice

3. DESTINATION OF USE

The OB pressure infusion cuff is a medical device and allows you to apply pressure (previously set) on the inserted infusion bag, this option is useful when, for example, the stand next to the patient (which uses gravity force) is not available or not possible to be used. This device is generally used during the moving of patients with stretchers, immobilization systems and specific transport solutions, even with pediatric stretchers.

The pressure infusion cuff has to be used only by medical or paramedical personnel, the field in which is applied is quite wide, from emergency situations to hospital or domestic use.

4. WARNING

The use of this device for purposes which are different from what stated in this document is discouraged. The Oscar Boscarol srl should not be considered responsible for damages caused by improper use or by a different use than what stated in this user manual. The manufacturer owns the right to change the device and the user manual without prior notice in order for functionality improvement.

Please refer to this manual for the correct use of the device.

NOTICE: Never exceed the pressure value of 300 mmHg: the device pressure directly affects the liquid infusion speed into the catheter. Such operation must be performed by trained and qualified personnel only.

5. ASSEMBLY

Connect rubber tubes from the cuff with the manometer and with the inflation bulb with regulating valve.

6. USE INSTRUCTIONS

Insert the infusion cuff into the built-in pouch, if any doubt arises refer to the picture available in the manual.

Close the valve situated on the bulb (clockwise) and inflate the cuff till the required pressure is reached (it is possible to control the applied pressure value on the manometer).

The value must be chosen depending on the medical substance which is to be injected and its required speed.

If necessary, release the faucet on the bulb to reach the required pressure.

In order to quickly deflate the cuff, open the faucet on the bulb entirely.

ATTENTION! Please always rely on qualified and trained personnel!

7. CLEANING

MANOMETER: Only use soft and dry cloth.

INFLATION BULB: Clean with soft and moist cloth and dry accurately (do NOT autoclave it)

FABRIC POUCH: Leave out the cuff from inside the pocket. Clean the cloth part with warm water and neutral soap (NO washing machine, NO spin)

BLADDER: Clean it with warm water and neutral soap paying attention that water does not permeate. Leave it to dry at room temperature

8. MAINTENANCE

Boscarol medical devices are accurately checked before commercialized and provided with CE mark.

Before every use, make sure the cuff is completely empty and verify that the pointer on the manometer display is at zero position with a tolerance of +/- mmHg.

Before utilize, in non-use conditions, always verify that pointer is acting inside tolerance area of ± 3 mmHg (recognizable by a small red line placed perpendicularly to the notches of graduate scale in proximity of zero value). Manometers out of ± 3 mmHg tolerance need to be recalibrated. Please refer to the manufacturer for technical assistance.

9. WARRANTY

The device is indicated for professional use and has a 12 months warranty from the date of purchase. The manufacturer, or the authorized service center, owns the right to repair or replace the parts affected by manufacturing or material faults, or to replace the entire device to its at their discretion.

Warranty is not applicable in case of damages caused by improper use, abuse, alteration or interventions performed by non-qualified unauthorized personnel.

With respect to this, please send the device to be repaired, together with proof of purchase, to the manufacturer or to the service assistance center stating the encounter failure/s.

The life of the device is 3 years if used and stored as indicated by the manufacturer.

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